



Laerdal
helping save lives



The Ideal Classroom Solution for CPR Training

Mini Anne Plus

Mini Anne Plus is a cost-effective solution for teaching quality CPR using cost effective durable and reusable manikins. Ten individual CPR manikins are included in each set, increasing student hands-on-time during training. The innovative new pump bag provides a simple, hygienic inflation method. Mini Anne Plus includes everything to get your CPR programme up and running in one convenient carry bag.

Mini Anne Plus



10 Mini Anne manikins, 10 kneel mats, 10 pump bags, 1 carry bag, 2 mesh collection bags, 10 face masks, 50 manikin wipes, 50 airways

FEATURES

- Reusable set of 10 Mini Anne manikins
- Durable, tested to at least 300,000 compressions (over 100 CPR classes)
- Easy-to-use pump bag for hygienic inflation
- Convenient storage and collection options
- Guidelines compliant clicker setting
- Soft clicker setting encourages new learners
- Carry bag, including pump bag and kneel mats for each student
- Quick preparation and clean up

ORDERING INFORMATION

Mini Anne Plus

- 106-00150 Mini Anne Plus 10 pack, light skin
106-00250 Mini Anne Plus 10 pack, dark skin

Includes: 10 manikins, 10 kneel mats, 10 pump bags, 1 carry bag, 2 mesh collection bags, 10 face masks, 50 manikin wipes, 50 airways

- 106-00350 Mini Anne Plus (single) light skin complete with Pump Bag
106-00450 Mini Anne Plus (single) dark skin complete with Pump Bag

Includes: 1 manikin, 1 kneel mat, 1 pump bag, 1 face mask, 2 manikin wipes, 1 airway

Consumables

- 106-10050 Airway Assembly MiniAnnePlus (pkg. 50)
106-11050 Face Mini Anne Light skin (pkg. 5)
106-12050 Face Mini Anne Dark skin (pkg. 5)
152400 Resusci Manikin Wipes
152401 Resusci Manikin Wipes

Accessories

- 106-10103 Instructor Pump (UK)
106-10101 Instructor Pump (US/JAP)
106-10115 Instructor Pump (Aus/CHN)
106-10133 Instructor Pump (Eur/KOR)
106-10250 Mesh Bag MiniAnnePlus (pkg. 2)
106-10350 Manikin MiniAnnePlus (single)Light skin
106-10450 Manikin MiniAnnePlus (single)Dark skin
106-10550 Kneel Mat (pkg. 10)
106-10650 Pump Bag MiniAnnePlus
103-00002 Mini Anne AED Training kit (pkg. 5)



Laerdal

helping save lives